NewSurg - Q & Pico

Q-Switched Laser Systems for the Safe And Effective Removal of Tattoos of All Colors, Pigmented Lesions and Dermal Toning

FDA Approved
High Rep Rate
Twin Pulse Mode
High Power Output
Multispot Zoom Handpiece
Pre-Programmed Parameters











NewSurg-Q

Specifications:

Wavelengths 1064 / 532 / 660 / 595

Pulse Duration 4 nanoseconds

Repetition Rate Up to 20 Hz

Spot Size Adjustable up to 10 mm

NewSurg-Pico

Specifications:

Wavelengths 1064 / 532 / 660 / 595

Pulse Duration 450 picoseconds

Repetition Rate Up to 10 Hz

Spot Size Adjustable up to 10 mm



Phone: (215) 822-7722 Email: sales@newsurg.com Website: NewSurg.com